

Tricefy and Athenahealth

Last Modified on 01/04/2021 6:00 am EST

Trice and athenahealth

Improving the way clinicians share with patients, collaborate with colleagues, and manage their image data

The seamless integration between Tricefy and athenaClinicals, has streamlined how users access their imaging data. Here's a short video to show you how it works:

1. Users place an order for an imaging study in athenaNet.
2. Next the study is performed and all pertinent DICOM data is sent to Tricefy.
3. Once completed, a link to the images will be available in the patient chart under "Imaging Results".
4. Clicking on the link labelled "trice images" will open the study into the Tricefy viewer.
5. All Tricefy viewer functions are available: sharing, collaboration, reporting

Your browser does not support HTML5 video.
